




Electro Surgical
General Surgical
Instruments
Catalog



ISO 9001 : 2015 - ISO 13485 : 2016, CE, FDA

www.asifvaincetrd.com



P.O. B.o.x. 3131
Commissioner Road Opposit Truama Center
Sialkot-51310 - Pakistan
Tel: +92-52-4566456
Mob: +92-333-8708000 / +92-331-6780029
E-mail: zahid@asifvaincetrd.com
Web: www.asifvaincetrd.com

Production Unit:
Ranger Road Near Kentt Housing
Langrewali, Sialkot-51310 - Pakistan

A heartfelt thank you to all our loyal customers, around the world. We are proud to present our latest surgical & electro surgical instruments cataloge. With this issue we have aggressively expanded into several new markets.

AVT is the forefront of technology for demonstrable precision Of quality and cost.

AVT has specialized in exactness of development Production and sales of electro & general surgical Instruments.

AVT dominates strength in the field of electro & general surgical instruments. All corporate offices; the teams of research & development, sales & marketing and of production & quality are in frequent interaction in family like environment to ensure that the company remains innovative and customers focused.

As usual, all our instruments are available for private label & custom-manufacturing. We thank you for our support and Look forward to man ears of continued success.

Sincerely ,

Muhammad Asif
(Chairman)



ASIF VAINCE TRADERS



ASIF VAINCE TRADERS

Electro Surgical Items



Asif Vaince Traders bipolar forceps are available in various lengths and shapes for different surgical and anatomical requirements. The reusable bipolar forceps make it easy to coagulate tissue. Sticking effects on the tips of the forceps are reduced to a minimum.

Advantages of Reusable European Bipolar Forceps

- *Highly polished tips reduce tissue adherence.*
- *Light closing force causes less fatigue over time.*
- *Nylon coating resists abrasion and provides superb insulation.*
- *Precision tip alignment aids fine control.*
- *Reusable, durable, and designed for multiple uses.*
- *The Instruments is autoclavable.*
- *Compatible with all common electrosurgery unit.*



AV-01

McPherson



AV-02

McPherson Ang.



AV-03

Jeweller #5 Str



AV-04

Jeweller #5 Ang



AV-05


Jeweller #7



AV-06

Adson

Tip Chart




Tip size can be changed as per your requirement, just add the following extensions
ft tip size at the end of Ref #: **TIP**- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-01	McPherson	0.5	10.2
AV-02	McPherson Angled	0.5	10.2
AV-03	Jeweller #5	0.5, 0.7	11.5
AV-04	Jeweller #5 Angled	0.5, 0.7	11.5
AV-05	Jeweller #7	0.5, 0.7	11.5
AV-06	Adson	1.0	12.0
AV-07	Adson	1.0	15.0




Tip Chart




Tip size can be changed as per your requirement, just add the following extensions
ft tip size at the end of Ref #: **TIP-** Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-011	Straight Angled	0.5, 1.0, 2.0	15.0
AV-012	Straight Angled	0.5, 1.0, 2.0	18.0
AV-013	Straight Angled	0.5, 1.0, 2.0	20.0
AV-014	Straight Angled	0.5, 1.0, 2.0	22.0
AV-015	Straight Angled	0.5, 1.0, 2.0	24.0




Tip Chart




Tip size can be changed as per your requirement, just add the following extensions
ft tip size at the end of Ref #: **TIP**- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-016	Bayonet	0.5, 1.0, 2.0	16.5
AV-017	Bayonet	0.5, 1.0, 2.0	18.0
AV-018	Bayonet	0.5, 1.0, 2.0	20.0
AV-019	Bayonet	0.5, 1.0, 2.0	22.0
AV-020	Bayonet	0.5, 1.0, 2.0	24.0




Tip Chart




Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: TIP- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-021	Bayonet Angled Up	0.5, 1.0, 2.0	16.5
AV-022	Bayonet Angled Up	0.5, 1.0, 2.0	18.0
AV-023	Bayonet Angled Up	0.5, 1.0, 2.0	20.0
AV-024	Bayonet Angled Up	0.5, 1.0, 2.0	22.0
AV-025	Bayonet Angled Up	0.5, 1.0, 2.0	24.0



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions
ft tip size at the end of Ref #: **TIP**- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm


Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-026	Bayonet Angled Down	0.5, 1.0, 2.0	16.5
AV-027	Bayonet Angled Down	0.5, 1.0, 2.0	18.0
AV-028	Bayonet Angled Down	0.5, 1.0, 2.0	20.0
AV-029	Bayonet Angled Down	0.5, 1.0, 2.0	22.0
AV-030	Bayonet Angled Down	0.5, 1.0, 2.0	24.0




It is well known that depending on the generator and the surface quality of the forceps jaws, coagulated tissue tends to more or less stick to the instrument jaw. This unwanted effect can occur particularly often with liver tissue and is critical in neurosurgery. For these applications in particular, we have designed our Bipolar Forceps with integrated irrigation channel. The coagulation area is irrigated through a built-in irrigation channel inside the forceps branch which can be connected to a standard irrigation hose. By means of this irrigation the sticking of the tissue is substantially reduced and surgical work is made easier.

Advantages of Reusable European Irrigation Bipolar Forceps

- *A universal hub connector allows irrigation tubing to be connected easily.*
- *The ergonomically designed handle provides less fatigue and a well balanced instrument.*
- *Nylon coating resists abrasion and provides superb insulation.*
- *For high visibility our irrigation tube can be recessed into the forceps blank (buried irrigation).*
- *Compatible with all common electrosurgery unit.*



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: TIP- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-031	Straight Irrigation	0.5, 1.0, 2.0	18.0
AV-032	Straight Irrigation	0.5, 1.0, 2.0	20.0
AV-033	Bayonet Irrigation	0.5, 1.0, 2.0	18.0
AV-034	Bayonet Irrigation	0.5, 1.0, 2.0	20.0

Note: More Bipolar Irrigation Forceps are available use IR extension to any regular bipolar forceps.




During the coagulation process, tissue and tips of forceps can stick together. **AVT** non-stick bipolar forceps are made with a special silver alloy and therefore prevent tissue from sticking. Disturbing and time consuming disruptions for cleaning the tips during surgery are thus avoided. Due to the good thermal qualities of silver the temperature of the forceps tips does not exceed 100°C during coagulation, which enhances the non-stick effect. With over ten years of experience we can guarantee that our bipolar forceps do not stick. The non-stick effect remains constant and is not reduced by frequent sterilization. Our non-stick bipolar forceps together with the appropriate connection cables are compatible with all common types of high-frequency units such as Erbe, Martin, Berchtold, Valleylab. **AVT** Bipolar forceps non-stick are available in straight, angled or bayonet shape. Available lengths range from 110 to 300 mm. Tips can be made from 0.25 to 2.0 mm width. Forceps are also available with different measurements as well as with flat pin or 2-pin plug.

Advantages of Reusable Non-Stick European Bipolar Forceps

- *Ten years of experience speak for themselves. Our forceps do not stick – never!*
- *No cleaning during surgery*
- *Especially high thermal qualities*
- *All Non-stick Bipolar Forceps are re-usable and autoclavable.*
- *Compatible with all common electrosurgery unit.*



Tip Chart




Tip size can be changed as per your requirement, just add the following extensions
ft tip size at the end of Ref #: TIP- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-035	McPherson	0.5,	10.2
AV-036	McPherson Angled	0.5,	10.2
AV-037	Jeweller #5	0.5, 0.7	11.5
AV-038	Jeweller #5 Angled	0.5, 0.7	11.5

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps.



Tip Chart




Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: TIP- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref		Tip Size (mm)	Length Size (cm)
AV-039	Straight	0.5, 1.0, 2.0	15.0
AV-040	Straight	0.5, 1.0, 2.0	18.0
AV-041	Straight	0.5, 1.0, 2.0	20.0
AV-042	Straight	0.5, 1.0, 2.0	22.0

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps.



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: TIP- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-043	Bayonet	0.5, 1.0, 2.0	16.5
AV-044	Bayonet	0.5, 1.0, 2.0	18.0
AV-045	Bayonet	0.5, 1.0, 2.0	20.0
AV-046	Bayonet	0.5, 1.0, 2.0	22.0


Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps.




Asif Vaince Traders bipolar forceps are available in various lengths and shapes for different surgical and anatomical requirements. The reusable Bipolar forceps make it easy to coagulate tissue. Sticking effects on the tips of the forceps are reduced to a minimum.

Advantages of Reusable US Bipolar Forceps

- *Highly polished tips reduce tissue adherence.*
- *Light closing force causes less fatigue over time.*
- *Nylon coating resists abrasion and provides superb insulation.*
- *Reusable, durable, and designed for multiple uses.*
- *Precision tip alignment aids fine control.*
- *The Instruments is atoclavable.*
- *Compatible with all common electrosurgery unit.*




Tip Chart




Tip size can be changed as per your requirement; just add the following extensions
At the end of the end of Ref. # 478: Needle, 0.5mm, 1.0mm, 1.5mm & 2.0mm.


Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-047	McPherson	0.5	10.2
AV-048	McPherson Angled	0.5	10.2
AV-049	Coaptation Nedler	0.5	10.2
AV-050	Coaptation	0.5	10.8
AV-051	Jeweller #5	Needle, 0.5	11.5




Tip Chart



Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-052	Jeweller #4	Needle, 0.5	11.5
AV-053	Tenzel	Needle, 0.5	10.2
AV-054	Jeweller #7	Needle, 0.5	11.5
AV-055	Iris Straight	Needle, 0.5	10.8
AV-056	Iris Angled	0.5, 1.0	10.8




Tip Chart




Tip size can be changed as per your requirement, just add the following extensions
At the end of the end of Ref #10 Needle: 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm


Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-057	Adson	1.0	12.0
AV-058	Semkin	0.5, 1.0	14.0
AV-059	Cushing	1.0, 1.5, 2.0	15.9
AV-060	Cushing	1.0, 1.5, 2.0	17.8



Tip Chart



Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-061	Cushing Gush	1.0, 1.5, 2.0	17.8
AV-062	Gerald Straight	1.0, 1.5, 2.0	17.8
AV-063	Gerald Angled	1.0, 1.5, 2.0	17.8
AV-064	Gerald Bayonet	1.0, 1.5, 2.0	18.4




Tip Chart


Needle	0.5mm	0.7mm	1.0mm	1.5mm	2.0mm

Tip size can be changed as per your requirement, just add the following extensions
to the cable at the end of Ref # - 000. Needle: 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-065	Jansen Bayonet	1.0, 1.5, 2.0	15.2
AV-066	Jansen Bayonet	1.0, 1.5, 2.0	17.8
AV-067	Jansen Bayonet	1.0, 1.5, 2.0	20.4
AV-068	Malis	1.0, 1.5, 2.0	17.8



Tip Chart



Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-069	Scoville-Greenwood	1.0, 1.5	19.1
AV-070	Yasargil	1.0, 1.5	19.7
AV-071	Hardy Bayonet	1.0, 1.5	21.6
AV-072	Hardy Bayonet Ang	1.0, 1.5	21.6



Our Monopolar instruments are tested at high volts to ensure proper functioning during the most stringent operations. With the best quality stainless steel, excellent finishing and superior insulation, we provide worldwide guarantee for all our Monopolar instruments.

Advantages of Reusable Monopolar Forceps

- *Lightweight and well balanced to prevent fatigue.*
- *Highly polished tips minimise sticking.*
- *Nylon coating resists abrasion and provides superb insulation.*
- *Precision tip alignment aids fine control.*
- *Compatible with all common electrosurgery unit.*

AV-073 Jeweller #4



AV-074 Jeweller #7



AV-075 Adson



AV-076 Potts-Smith



AV-077 Potts-Smith



AV-078 Jansen Bayonet




AV-079 Jansen Bayonet



AV-080 Jansen Bayonet UP



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions
ft tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-073	Jeweller #4	Needle	11.5
AV-074	Jeweller #7	Needle	11.5
AV-075	Adson	1.0	12.0
AV-076	Potts-Smith	2.0	17.8
AV-077	Potts-Smith	2.0	21.0
AV-077-1	Potts-Smith	2.0	24.1
AV-078	Jansen Bayonet	1.5	13.3
AV-078-1	Jansen Bayonet	1.5	15.9
AV-079	Jansen Bayonet	2.0	19.1
AV-079-1	Jansen Bayonet	2.0	21.0
AV-080	Jansen Bayonet UP	2.0	19.1
AV-080-1	Jansen Bayonet UP	2.0	21.0



It is well known that depending on the generator and the surface quality of the forceps jaws, coagulated tissue tends to more or less stick to the instrument jaw. This unwanted effect can occur particularly often with liver tissue and is critical in neurosurgery. For these applications in particular, we have designed our Bipolar Forceps with integrated irrigation channel. The coagulation area is irrigated through a built-in irrigation channel inside the forceps branch which can be connected to a standard irrigation hose. By means of this irrigation the sticking of the tissue is substantially reduced and surgical work is made easier.

Advantages of Reusable US Irrigation Forceps

- *A universal hub connector allows irrigation tubing to be connected easily.*
- *The ergonomically designed handle provides less fatigue and a well balanced instrument.*
- *For high visibility our irrigation tube can be recessed into the forceps blank (buried irrigation).*
- *Compatible with all common electrosurgery unit.*



AV-o81

Hardy Bayonet



AV-o82

Jansen Bayonet



AV-o83

Scoville-Greenwood



AV-o84


Gerald Bayonet



AV-o85

Straight

Tip Chart



Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-o81 -IR	Hardy Bayonet	0.5, 1.0, 1.5	21.6
AV-o82 -IR	Jansen Bayonet	0.5, 1.0, 1.5	19.1
AV-o83 -IR	Scoville-Greenwood	0.5, 1.0, 1.5	19.7
AV-o84 -IR	Gerald Bayonet	0.5, 1.0, 1.5	18.4
AV-o85 -IR	Straight	0.5, 1.0, 1.5	15.9

Note: More Irrigation Bipolar Forceps are available use IR extension to any regular bipolar forceps.



During the coagulation process, tissue and tips of forceps can stick together. **Asif** non-stick bipolar forceps are made with a special silver alloy and therefore prevent tissue from sticking. Disturbing and time consuming disruptions for cleaning the tips during surgery are thus avoided. Due to the good thermal qualities of silver the temperature of the forceps tips does not exceed 100°C during coagulation, which enhances the non-stick effect. With over ten years of experience we can guarantee that our bipolar forceps do not stick. The non-stick effect remains constant and is not reduced by frequent sterilization. Our non-stick bipolar forceps together with the appropriate connection cables are compatible with all common types of high-frequency units such as Erbe, Martin, Berchtold, Valleylab. **Asif** Bipolar forceps non-stick are available in straight, angled or bayonet shape. Available lengths range from 110 to 300 mm. Tips can be made from 0.25 to 2.0 mm width. Forceps are also available with different measurements as well as with flat pin or 2-pin plug.

Advantages of Reusable Non-Stick European Bipolar Forceps

- *Ten years of experience speak for themselves. Our forceps do not stick – never!*
- *No cleaning during surgery*
- *Especially high thermal qualities*
- *All Non-stick Bipolar Forceps are re-usable and autoclavable.*
- *Compatible with all common electrosurgery unit.*



AV-086

McPherson



AV-087

McPherson Angled



AV-088

Jeweller #5



AV-089

Jeweller #4



AV-090


Adson



AV-091

Semkin

Tip Chart



Tip size can be changed as per your requirement, just add the following extensions
for tip size at the end of Ref #: TIP- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-086	McPherson	0.5	10.2
AV-087	McPherson Angled	0.5	10.2
AV-088	Jeweller #5	Needle, 0.5	11.5
AV-089	Jeweller #4	Needle, 0.5	11.5
AV-090	Adson	1.0	12.0
AV-091	Semkin	0.5, 1.0	14.0

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps



AV-092

Cushing



AV-093

Gerald Bayonet



AV-094

Jansen Bayonet



AV-095

Scoville-Greenwood



AV-096


Yasargil



AV-097

Hardy Bayonet

Tip Chart



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV - 092	Cushing	1.0, 1.5, 2.0	15.9
AV - 092	Cushing	1.0, 1.5, 2.0	17.8
AV - 092	Gerald Bayonet	1.0, 1.5	18.4
AV - 092	Jansen Bayonet	1.0, 1.5	15.9
AV - 092	Scoville-Greenwood	1.0, 1.5	19.1
AV - 092	Yasargil	1.0, 1.5	19.7
AV - 092	Hardy Bayonet	1.0, 1.5	21.6

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps



AV-092

Cushing



AV-093

Gerald Bayonet



AV-094

Jansen Bayonet



AV-095

Scoville-Greenwood



AV-096


Yasargil



AV-097

Hardy Bayonet



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: **TIP**- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV - 092	Cushing	1.0, 1.5, 2.0	15.9
AV - 092	Cushing	1.0, 1.5, 2.0	17.8
AV - 092	Gerald Bayonet	1.0, 1.5	18.4
AV - 092	Jansen Bayonet	1.0, 1.5	15.9
AV - 092	Scoville-Greenwood	1.0, 1.5	19.1
AV - 092	Yasargil	1.0, 1.5	19.7
AV - 092	Hardy Bayonet	1.0, 1.5	21.6

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps



Tip size can be changed as per your requirement, just add the following extensions
 fit tip size at the end of Ref #: **TIP-** Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-098	Jeweller #4	0.5,	11.5
AV-099	Semkin	0.7	14.0
AV-100	Jansen Bayonet	0.7	19.1
AV-101	Scoville-Greenwood	1.5	19.7



AV-102

Jansen Bayonet



AV-103

Jansen Bayonet



AV-104

Potts-Smith




AV-105

Potts-Smith



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-102	Jansen Bayonet	2.0	19.0
AV-103	Jansen Bayonet	2.0	21.0
AV-104	Potts-Smith	2.0	21.0
AV-105	Potts-Smith	2.0	24.0



Our Diathermy instruments are tested at high volts to ensure proper functioning during the most stringent operations. With the best quality stainless steel, excellent finishing and superior insulation, we provide worldwide guarantee for all our Diathermy instruments. Diathermy instruments are available in four base types. Models popular Diathermy Instruments : Adson, Jeweller, McIndoe, Standard, Bayonet, Artery Forceps, Wilson Hey Forceps, Diathermy Scissors, etc in various lengths, tip sizes, and 4 plug types

Advantages of Reusable Diathermy Instruments

- *Lightweight and well balanced to prevent fatigue.*
- *Highly polished tips minimise sticking.*
- *Precision tip alignment aids fine control.*
- *Compatible with all common electrosurgery unit.*



Tuner Warwick Scissors (STR)
AV-106



Tuner Warwick Scissors (CRV)
AV-107



Dressing Scissors (STR)
AV-110



Dressing Scissors (CRV)
AV-111



Wilson Hey Forceps (STR)
AV-114



Wilson Hey Forceps (CRV)
AV-115



Artery Forceps (STR)
AV-118



Artery Forceps (CRV)
AV-119



Abbey Needle



Richie Forceps
AV-138

ENGLISH MONOPOLAR MALE 4.8MM CONNECTOR

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-106	Tuner Warwick Scissors (STR)	-	17.8
AV-107	Tuner Warwick Scissors (CRV)	-	17.8
AV-108	Tuner Warwick Scissors (STR)	-	23.0
AV-109	Tuner Warwick Scissors (CRV)	-	23.0
AV-110	Dressing Scissors (STR)	-	17.8
AV-111	Dressing Scissors (CRV)	-	17.8
AV-112	Dressing Scissors (STR)	-	23.0
AV-113	Dressing Scissors (CRV)	-	23.0
AV-114	Wilson Hay Forceps (STR)	-	17.8
AV-115	Wilson Hay Forceps (CRV)	-	17.8
AV-116	Wilson Hay Forceps (STR)	-	23.0
AV-117	Wilson Hay Forceps (CRV)	-	23.0
AV-118	Artery Forceps (STR)	-	17.8
AV-119	Artery Forceps (CRV)	-	17.8
AV-120	Artery Forceps (STR)	-	23.0
AV-121	Artery Forceps (CRV)	-	23.0

EUROPEAN MONOPOLAR FEMALE 4.0MM CONNECTOR

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-122	Tuner Warwick Scissors (STR)	-	17.8
AV-123	Tuner Warwick Scissors (CRV)	-	17.8
AV-124	Tuner Warwick Scissors (STR)	-	23.0
AV-125	Tuner Warwick Scissors (CRV)	-	23.0
AV-126	Dressing Scissors (STR)	-	17.8
AV-127	Dressing Scissors (CRV)	-	17.8
AV-128	Dressing Scissors (STR)	-	23.0
AV-129	Dressing Scissors (CRV)	-	23.0
AV-130	Wilson Hay Forceps (STR)	-	17.8
AV-131	Wilson Hay Forceps (CRV)	-	17.8
AV-132	Wilson Hay Forceps (STR)	-	23.0
AV-133	Wilson Hay Forceps (CRV)	-	23.0
AV-134	Artery Forceps (STR)	-	17.8
AV-135	Artery Forceps (CRV)	-	17.8
AV-136	Artery Forceps (STR)	-	23.0
AV-137	Artery Forceps (CRV)	-	23.0



Adson
AV-139



McINDOE
AV-140



Gillies
AV-141



Dissecting
AV-145



Gerald
AV-149



Adson Bayonet
AV-150



Adson Bayonet UP
AV-152



Remington Hobb
AV-155

ENGLISH MONOPOLAR MALE 4.8MM CONNECTOR

Ref	Item Name	Tip Size	Length Size
AV-139	Adson	1.0	12.0
AV-140	McINDOE	1.0, 2.0	15.0
AV-141	Gillies	-	15.0
AV-142	Gillies	-	17.8
AV-143	Gillies	-	20.0
AV-144	Gillies	-	24.0
AV-145	Dissecting	1.0, 2.0	15.0
AV-146	Dissecting	1.0, 2.0	17.8
AV-147	Dissecting	1.0, 2.0	21.0
AV-148	Dissecting	1.0, 2.0	24.1
AV-149	Gerald 1x2	-	17.8
AV-150	Adson Bayonet	1.0, 2.0	16.0
AV-151	Adson Bayonet	1.0, 2.0	19.1
AV-152	Adson Bayonet Angled UP	1.0, 2.0	16.0
AV-153	Adson Bayonet Angled UP	1.0, 2.0	19.1
AV-154	Adson Bayonet Angled UP	1.0, 2.0	21.0
AV-155	Remington Hobb	1.0, 2.0	24.1

EUROPEAN MONOPOLAR FEMALE 4.0MM CONNECTOR

Ref	Item Name	Tip Size	Length Size
AV-156	Adson	1.0	12.0
AV-157	McINDOE	1.0, 2.0	15.0
AV-158	Gillies	-	15.0
AV-159	Gillies	-	17.8
AV-160	Gillies	-	20.0
AV-161	Gillies	-	24.0
AV-162	Dissecting	1.0, 2.0	15.0
AV-163	Dissecting	1.0, 2.0	17.8
AV-164	Dissecting	1.0, 2.0	21.0
AV-165	Dissecting	1.0, 2.0	24.1
AV-166	Gerald 1x2	-	17.8
AV-167	Adson Bayonet	1.0, 2.0	16.0
AV-168	Adson Bayonet	1.0, 2.0	19.1
AV-169	Adson Bayonet Angled UP	1.0, 2.0	16.0
AV-170	Adson Bayonet Angled UP	1.0, 2.0	19.1
AV-171	Adson Bayonet Angled UP	1.0, 2.0	21.0
AV-172	Remington Hobb	1.0, 2.0	24.1

ENGLISH MONOPOLAR MALE 4.0MM CONNECTOR

Ref	Item Name	Tip Size	Length Size
AV-173	Adson	1.0	12.0
AV-174	McINDOE	1.0, 2.0	15.0
AV-175	Gillies	-	15.0
AV-176	Gillies	-	17.8
AV-177	Gillies	-	20.0
AV-178	Gillies	-	24.0
AV-179	Dissecting	1.0, 2.0	15.0
AV-180	Dissecting	1.0, 2.0	17.8
AV-181	Dissecting	1.0, 2.0	21.0
AV-182	Dissecting	1.0, 2.0	24.1
AV-183	Gerald 1x2	-	17.8
AV-184	Adson Bayonet	1.0, 2.0	16.0
AV-185	Adson Bayonet	1.0, 2.0	19.1
AV-186	Adson Bayonet Angled UP	1.0, 2.0	16.0
AV-187	Adson Bayonet Angled UP	1.0, 2.0	19.1
AV-188	Adson Bayonet Angled UP	1.0, 2.0	21.0
AV-189	Remington Hobb	1.0, 2.0	24.1

MONOPOLAR WITHOUT CONNECTOR

Ref	Item Name	Tip Size	Length Size
AV-190	Adson	1.0	12.0
AV-191	McINDOE	1.0, 2.0	15.0
AV-192	Gillies	-	15.0
AV-193	Gillies	-	17.8
AV-194	Gillies	-	20.0
AV-195	Gillies	-	24.0
AV-196	Dissecting	1.0, 2.0	15.0
AV-197	Dissecting	1.0, 2.0	17.8
AV-198	Dissecting	1.0, 2.0	21.0
AV-199	Dissecting	1.0, 2.0	24.1
AV-200	Gerald 1x2	-	17.8
AV-201	Adson Bayonet	1.0, 2.0	16.0
AV-202	Adson Bayonet	1.0, 2.0	19.1
AV-203	Adson Bayonet Angled UP	1.0, 2.0	16.0
AV-204	Adson Bayonet Angled UP	1.0, 2.0	19.1
AV-205	Adson Bayonet Angled UP	1.0, 2.0	21.0
AV-206	Remington Hobb	1.0, 2.0	24.1



Asif Vaince Traders disposable bipolar forceps have been designed and tested to meet the practical. The main purpose for offering single-use bipolar forceps is to avoid the risk of contamination while remaining cost effective and utilizing the functionalities of the higher end forceps. Hardened stainless steel provides an optimum spring tension and serves as the basis for a superior instrument.

Advantages of Disposable Bipolar Forceps

- *Lightweight and well balanced to prevent fatigue.*
- *Highly polished tips minimise sticking.*
- *Nylon coating resists abrasion and provides superb insulation.*
- *Precision tip alignment aids fine control.*
- *Cost effective in comparison to re-usable.*
- *All Bipolar forceps available with cable, sterile packed.*
- *Compatible with all common electrosurgery unit.*



AV-207

Jeweller #5



AV-208

Straight



AV-209

Bayonet



AV-210


Straight



AV-211

Straight Angled

Tip Chart




Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: TIP- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-207	Jeweller #5 With 3m Cable	0.5	11.5
AV-208	Straight With 3m Cable	1.0, 1.5	18.0
AV-209	Bayonet With 3m Cable	1.0, 1.5	18.0
AV-210	Straight	1.0, 1.5	18.0
AV-211	Bayonet	1.0, 1.5	20.0

Note: More Disposable Bipolar Forceps are Available Use DP extension to any regular Bipolar Forceps.



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions
for tip size at the end of Ref #: **TIP**- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
AV-212	US. Gerald	1.0, 1.5	17.8
AV-213	US. Jansen Bayonet Irrigation	1.0, 1.5	19.1
AV-214	US. Jansen Bayonet	1.0, 1.5	17.8
AV-215	US. Jeweller #5	0.5	11.5
AV-216	US. Adson	1.0	12.1

Note: More Disposable Bipolar Forceps are Available Use DP extension to any regular Bipolar Forceps.



The First Choice For Proven Performance, **Asif Vaince Traders** Provides The Broadest selection of fully autoclavable instruments for LEEP procedures. Our "rose coated" instruments will provide you with durability and performance unmatched in the industry. Models Available: Graves, Pederson, VU-More, 4-Way, Cusco, Collin, Weisman, Gelpi, Rigby, Lateral, Lateral Vaginal, Endocervical, Iris, Emmet, Straight, two-prong, three-prong, etc.



Graves Vaginal Speculum

With Smoke Evacuation Tube

Small AV-217

Medium AV-218

Large AV-219

Disposable Tube



Small AV-220

Medium AV-221

Large AV-222



Pederson Vaginal Speculum

With Smoke Evacuation Tube

Small AV-223

Medium AV-224

Large AV-225

Disposable Tube



Small AV-226

Medium AV-227

Large AV-228



VU-Max Grave Vaginal Speculum

With Smoke Evacuation Tube

Medium AV-229

Large AV-230

Disposable Tube



Medium AV-231

Large AV-232



VU-Max Pederson Vaginal Speculum

With Smoke Evacuation Tube

Medium AV-233

Large AV-234

Disposable Tube



Medium AV-235

Large AV-236



Collin Vaginal Speculum

With Smoke Evacuation Tube

Small AV-237

Medium AV-238



CUSCO Vaginal Speculum

With Smoke Evacuation Tube

Small AV-240

Medium AV-241



Disposable Tube



VU-More Graves Vaginal Speculum
With Smoke Evacuation Tube
Medium AV-245
Large AV-246



4-Blades Speculum
With Smoke Evacuation Tube
Medium AV-247
Large AV-248



Disposable Tube



**Weisman Vaginal Speculum
Right Side Open**
With Smoke Evacuation Tube
Medium AV-249
Large AV-250



Disposable Tube



**Weisman Vaginal Speculum
Left Side Open**
With Smoke Evacuation Tube
Medium AV-253
Large AV-254

Medium AV-255
Large AV-256



**Lateral Retractor
Wide Open**
AV-257



**Rigby Retractor
Wide Open**
AV-258



**Gelpi Retractor
Wide Open**
AV-259



Kogan Speculum
3mm AV-260
4mm AV-261
5mm AV-262



**Lateral Retractor
21cm 22mmx75mm**
AV-263



Schroeder
Tanaculum 24cm

AV-264



Schroeder
Tanaculum 24cm

AV-265



Single Tooth
Tanaculum 24cm

AV-267



Campion
Forceps 22cm

AV-268



Spong
Forceps 22cm

AV-269

Dressing Forceps

AV-270 21cm

AV-271 25cm



Tissue Forceps

AV-272 21cm

AV-273 25cm



Iris Hook

AV-274 25cm



Emmet Hook

AV-275 25cm



Emmet Hook Straight

AV-276 25cm



Emmet Hook Two Prong

AV-277 25cm



Emmet Hook Three Prong

AV-278 25cm





Through a non-reflective coating technique the glare produced during laser surgery is significantly minimized. The stainless steel surface is specially prepared keeping in mind the superior coating requirements for our range of laser surgery instruments. Custom sizes and shapes are available as well as partial gold plating services for certain instruments as per customer requirements. Models Available: Graves, Pederson, VU-More, 4-Way, Cusco, Collin, Weisman, Gelpi, Rigby, Lateral, Lateral Vaginal, Endocervical, Iris, Emmet, Straight, two-prong, three-prong, etc.



Graves Vaginal Speculum

With Smoke Evacuation Tube

Small **AV-279**

Medium **AV-280**

Large **AV-281**



Pederson Vaginal Speculum

With Smoke Evacuation Tube

Small **AV-283**

Medium **AV-284**

Large **AV-285**




Collin Vaginal Speculum

With Smoke Evacuation Tube

Small **AV-287**

Medium **AV-288**

Large **AV-289**



Lateral Vaginal Speculum

AV-282



Kogan Endocervical Speculum

AV-286



3mm 381-770



4mm 381-772



5mm 381-774



**Rigby Retractor
AV-290**



**Gelpi Retractor
AV-291**



**Schroeder Tenaculum
AV-292**



**Iris Hook
AV-293**



**Emmet Hook
AV-294**



**Emmet Hook Straight
AV-295**



**Emmet Hook Two Prong
AV-296**





These E-Z Clean Laparoscopic Electrodes are made from very high quality raw material which ensures hassle free work performance at its user end. These E-Z Clean Laparoscopic Electrodes are widely finds its applications in various medical sectors.

Salient Features of our E-Z Clean Laparoscopic Electrodes

- *Provide excellent dissection and coagulation for a variety of surgical applications.*
- *High-tech PTFE coating for reduced eschar buildup.*
- *Reintroduction during surgery is minimized.*
- *Lower power settings for safer delivery of electrosurgery.*
- *Tips are cost-effective because they Save Time and Money.*



Without Suction



L-Hook Electrode 36cm
AV-298



J-Hook Electrode 36cm
AV-299



Needle Electrode 36cm
AV-300



Spectula Electrode 36cm
AV-301



Ball Electrode 36cm
AV-302



With Suction



L-Hook Electrode 36cm
AV-303



J-Hook Electrode 36cm
AV-304



Needle Electrode 36cm
AV-305



Spectula Electrode 36cm
AV-306



Ball Electrode 36cm
AV-307



Dagger 12xm
AV-308



Dagger 18cm
AV-309



Hook 12cm
AV-310



Hook 18cm
AV-311



Hockey 12cm
AV-312



Hockey 18cm
AV-313



Frazier Suction Tube 11FR 13.5cm
AV-314



Frazier Suction Tube 11FR 20.5cm
AV-315



Frazier Suction Tube 11FR 40.6cm
AV-316





Electrode handle with blade electrode for single-use. Sterile packaging in individual pouches. Ready to use. This handle is made of high quality material and fits into the surgeon's hand comfortably due to its ergonomic shape. The quality meets highest expectations and fulfills all required international safety regulations.

- Single-use electrode handle including blade electrode and 3 m cable (also available: 5 m)
- Individually packaged, sterile.
- Push buttons for cutting and coagulation.
- Unit side: 3-pin connector.
- Tips are cost-effective because they Save Time and Money.
- Suitable also for use with ERBE, Martin or Berchtold units, by addition of an adapter especially designed for such purpose.



Rocker Switch Pencil (Single Use)
3m Cable AV-317



Foot Control Pencil (Single Use)
3m Cable AV-318
5m Cable AV-319



Finger Switch Pencil (Reusable)
3m Silicon Cable AV-320



Electrode Handle (Reusable)




Finger Switch Pencil (Single Use)
3m Cable AV-322
5m Cable AV-323



Blade Electrode

Disposable: AV-324-s Reusable: AV-324-R



Blade Electrode Crv

Disposable: AV-325-s Reusable: AV-325-R



Needle Electrode

Disposable: AV-326-s Reusable: AV-326-R



Needle Electrode Crv

Disposable: AV-327-s Reusable: AV-327-R



Ball Electrode 2mm

Disposable: AV-328-s Reusable: AV-328-R



Ball Electrode Crv 2mm

Disposable: AV-329-s Reusable: AV-329-R



Ball Electrode 3mm

Disposable: AV-330-s Reusable: AV-330-R



Ball Electrode 3mm Crv

Disposable: AV-331-s Reusable: AV-331-R



Ball Electrode 4mm

Disposable: AV-332-s Reusable: AV-332-R



Ball Electrode 4mm Crv

Disposable: AV-333-s Reusable: AV-333-R



Ball Electrode 5mm

Disposable: AV-343-s Reusable: AV-334-R



Ball Electrode 5mm Crv



Blade Electrode Shaft 13cm

Disposable: AV-336-s Reusable: AV-336-R



Needle Electrode Shaft 13cm

Disposable: AV-337-s Reusable: AV-337-R



Ball Electrode 3mm Shaft 13cm

Disposable: AV-338-s Reusable: AV-338-R



Ball Electrode 5mm Shaft 13cm

Disposable: AV-339-s Reusable: AV-339-R



Conization Electrode

Disposable: AV-340-s Reusable: AV-340-R



Circle Loop Electrode 4mm

Disposable: AV-341-s Reusable: AV-341-R



Circle Loop Electrode 6mm

Disposable: AV-342-s Reusable: AV-342-R



Circle Loop Electrode 8mm

Disposable: AV-343-s Reusable: AV-343-R



Diamond Electrode






Disposable: AV-344-s Reusable: AV-344-R








Fine Needle Electrode

Disposable: AV-345-s Reusable: AV-345-R




					
Size	1.0x0.4cm	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm
Disposable	AV-346-S	AV-347-S	AV-348-S	AV-349-S	AV-350-S
Reusable	AV-346-R	AV-347-R	AV-348-R	AV-349-R	AV-350-R









					
Size	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm	1.0x1.5cm
Disposable	AV-351-S	AV-352-S	AV-353-S	AV-354-S	AV-355-S
Reusable	AV-351-R	AV-352-R	AV-353-R	AV-354-R	AV-355-R



						
Size	1.5x0.5cm	1.5x0.7cm	1.5x0.8cm	1.5x1.0cm	1.5x1.2cm	1.5x1.5cm
Disposable	AV-356-S	AV-357-S	AV-358-S	AV-359-S	AV-360-S	AV-361-S
Reusable	AV-356-R	AV-357-R	AV-358-R	AV-359-R	AV-360-R	AV-361-R



						
Size	2.0x0.8cm	2.0x1.0cm	2.0x1.2cm	2.0x1.3cm	2.0x1.5cm	2.0x2.0cm
Disposable	AV-362-S	AV-363-S	AV-364-S	AV-365-S	AV-367-S	AV-368-S
Reusable	AV-362-R	AV-363-R	AV-364-R	AV-365-R	AV-367-R	AV-368-R



Blade Electrode

Disposable :AV-369-S Reusable:AV-369-R



Blade Electrode Crv

Disposable :AV-370-S Reusable:AV-370-R



Needle Electrode

Disposable :AV-371-S Reusable:AV-371-R



Needle Electrode Crv

Disposable :AV-372-S Reusable:AV-372-R



Ball Electrode 2mm

Disposable :AV-373-S Reusable:AV-373-R



Ball Electrode Crv 2mm

Disposable :AV-374-S Reusable:AV-374-R



Ball Electrode 3mm

Disposable :AV-375-S Reusable:AV-375-R



Ball Electrode 3mm Crv

Disposable :AV-376-S Reusable:AV-376-R



Ball Electrode 4mm

Disposable :AV-377-S Reusable:AV-377-R



Ball Electrode 4mm Crv

Disposable :AV-378-S Reusable:AV-378-R



Ball Electrode 5mm

Disposable :AV-389-S Reusable:AV-389-R



Ball Electrode 5mm Crv

Disposable :AV-390-S Reusable:AV-390-R



Blade Electrode Shaft 13cm

Disposable :AV-391-S Reusable:AV-391-R



Needle Electrode Shaft 13cm

Disposable :AV-392-S Reusable:AV-392-R



Ball Electrode Shaft 3mm 13cm

Disposable :AV-393-S Reusable:AV-393-R



Ball Electrode 4mm Shaft 13cm

Disposable :AV-394-S Reusable:AV-394-R



Ball Electrode 5mm 13cm

Disposable :AV-395-S Reusable:AV-395-R



Circle Loop Electrode 4mm 0.7cm

Disposable :AV-396-S Reusable:AV-396-R



Circle Loop Electrode 6mm 0.7cm

Disposable :AV-396-S Reusable:AV-396-R



Circle Loop Electrode 8mm 0.7cm

Disposable :AV-397-S Reusable:AV-397-R



Circle Loop Electrode 10mm 0.7cm

Disposable :AV-398-S Reusable:AV-398-R



Plate Electrode Small






Disposable :AV-399-S Reusable:AV-399-R








Plate Electrode Large

Disposable :AV-400-S Reusable:AV-400-R









					
Size	1.0x0.4cm	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm
Disposable	AV-401-S	AV-402-S	AV-403-S	AV-404-S	AV-405-S
Reusable	AV-401-R	AV-402-R	AV-403-R	AV-404-R	AV-405-R









					
Size	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm	1.0x1.5cm
Disposable	AV-406-S	AV-407-S	AV-408-S	AV-409-S	AV-410-S
Reusable	AV-406-R	AV-407-R	AV-408-R	AV-409-R	AV-410-R



						
Size	1.5x0.5cm	1.5x0.7cm	1.5x0.8cm	1.5x1.0cm	1.5x1.2cm	1.5x1.5cm
Disposable	AV-411-S	AV-412-S	AV-413-S	AV-414-S	AV-415-S	AV-416-S
Reusable	AV-411-R	AV-412-R	AV-413-R	AV-414-R	AV-415-R	AV-416-R



						
Size	2.0x0.8cm	2.0x1.0cm	2.0x1.2cm	2.0x1.3cm	2.0x1.5cm	2.0x2.0cm
Disposable	AV-417-S	AV-418-S	AV-419-S	AV-420-S	AV-421-S	AV-422-S
Reusable	AV-417-R	AV-418-R	AV-419-R	AV-420-R	AV-421-R	AV-422-R



Asif Vaince Traders Cables are available in disposable and reusable grade. We offer high quality silicon cables even for our single-use instruments thus making sure that all cables are able to withstand electrical current above the standard requirements. Furthermore, numerous tests including continuity and leak test ensures that our cables will work under the toughest requirements. Multiples sizes and plugs are available as per the need of the customer. Our cables are durable, flexible and able to withstand sterilization processes if needed. By following international standards we ensure that our cables are compatible with most major ESU's. Our cables are available in lengths 3 m, 4 m or 5 m.



Bipolar Cable

Instruments Connector: USA 2-Pin Plug

Generator Connector: 4mm Banana Plug

Reusable Silicon 3m Cable **Av-423**

Single Use 3m Cable **Av-424**



Bipolar Cable

Instruments Connector: European Flat Plug

Generator Connector: 4mm Banana Plug

Reusable Silicon 3m Cable **Av-425**

Single Use 3m Cable **Av-426**



Bipolar Cable

Instruments Connector: USA 2-Pin Plug

Generator Connector: Pin Distance 29mm

Reusable Silicon 3m Cable **Av-427**

Single Use 3m Cable **Av-428**



Bipolar Cable

Instruments Connector: European Flat Plug

Generator Connector: Pin Distance 29mm

Reusable Silicon 3m Cable **Av-429**

Single Use 3m Cable **Av-430**



Bipolar Cable

Instruments Connector: Bipolar Artery Sealer

Generator Connector: 4mm Banana Plug

Reusable Silicon 3m Cable **Av-431**



European Bipolar Cable

Instruments Connector: **European Flat Plug**

Generator Connector: **Martin Type**

Reusable Silicon 3m Cable **Av-432**



US Bipolar Cable

Instruments Connector: **USA 2-Pin Plug**

Generator Connector: **Martin Type**

Reusable Silicon 3m Cable **Av-433**



European Bipolar Cable

Instruments Connector: **European Flat Plug**

Generator Connector: **Erbe Type**

Reusable Silicon 3m Cable **Av-434**



US Bipolar Cable

Instruments Connector: **USA 2-Pin Plug**

Generator Connector: **Erbe Type**

Reusable Silicon 3m Cable **Av-435**



Bipolar Cable

Instruments Connector: **Ethicon**

Generator Connector: **4mm Banana Plug**

Reusable Silicon 3m Cable **Av-436**



Monopolar Cable

Instruments Connector: **Screw Type**

Generator Connector: **Standard Bovie 8mm Plug**

Reusable Silicon 3m Cable **AV-437**



Monopolar Cable

Instruments Connector: **4mm Female Fitting**

Generator Connector: **Standard Bovie 8mm Plug**

Reusable Silicon 3m Cable **AV-438**



Monopolar Endoscopic Cable

Instruments Connector: **4mm Female Fitting**

Generator Connector: **4mm Banana Plug**

Reusable Silicon 3m Cable **AV-439**



Diathermy Cable

Instruments Connector: **4mm Male Fitting**

Generator Connector: **4mm Banana Plug**

Reusable Silicon 3m Cable **AV-440**



Diathermy Cable

Instruments Connector: **4.8mm Female Fitting**

Generator Connector: **4mm Banana Plug**

Reusable Silicon 3m Cable **AV-441**



Diathermy Cable

Instruments Connector: **4.8mm Female Fitting**

Generator Connector: **Stan. Bovie 8mm Plug**

Reusable Silicon 3m Cable **AV-442**



- AV-443** Artery Sealer 6.5"
- AV-444** Artery Sealer 7.5"
- AV-445** Artery Sealer 8.5"
- AV-446** Artery Sealer 9.5"
- AV-447** Artery Sealer 10.5"



- AV-448** Non-Stick Artery Sealer 6.5"
- AV-449** Non-Stick Artery Sealer 7.5"
- AV-450** Non-Stick Artery Sealer 8.5"
- AV-451** Non-Stick Artery Sealer 9.5"
- AV-452** Non-Stick Artery Sealer 10.5"





ASIF VAINCE TRADERS

General Surgical Items



MAYO SCISSOR

MAYO SCISSOR

MAYO DISSECTING SCISSORS - STRAIGHT 6-3/4"
OVERALL LENGTH 6-3/4" (171.5MM)
STANDARD GRADE STAINLESS STEEL
MAYO DISSECTING SCISSORS
STRAIGHT END WITH SLIGHTLY BLUNT TIPS

METZENBAUM SCISSOR

ESSENTIAL SCISSORS FOR CUTTING DELICATE TISSUE AND BLUNT DISSECTION. HAVING A LONGER HANDLE-TO-BLADE DISTANCE, THE METZENBAUM IS AVAILABLE IN A VARIETY OF LENGTHS; THE 7-INCH LENGTH IS MOST PREFERRED.

METZENBAUM SCISSOR



IRIS SCISSOR

IRIS SCISSOR 4.5

REMOVE NECROTIC TISSUE
AVAILABLE IN STRAIGHT OR CURVED STYLES
CURVED SCISSORS FOLLOW TISSUE CONTOURS IN SMALL AREAS; STRAIGHT SCISSORS CUT ALONG FLAT SURFACES

DEAN SCISSOR

DEAN STANDARD SERIES SCISSORS WITH 1 SERRATED BLADE USED FOR TRIMMING TISSUE OR CUTTING SUTURES. 16.5CM



DEAN SCISSOR

Mayo- Heger Needle holder




MEDGYN CARRIES A COMPLETE LINE OF HIGH QUALITY NEEDLE HOLDERS.

ADDITIONAL PRODUCT INFORMATION

MAYO-HEGAR NEEDLE HOLDER
127MM SERRATED JAWS
152MM SERRATED JAWS
178MM SERRATED JAWS
203MM SERRATED JAWS

MOSQUITO FORCEPS 5"
HEMOSTATIC FORCEPS,
HALSTEAD MOSQUITO - CURVED 5"
> OVERALL LENGTH 5" (127MM)
> STANDARD GRADE STAINLESS STEEL
HALSTEAD MOSQUITO FORCEPS
> ADJUSTABLE LOCKING CROSS ARMS FOR
TIGHT CLAMPING, FINGER RING GNP
> CURVED HEMOSTATIC FORCEPS WITH
SERRATED TIPS



MOSQUITO FORCEPS



SPOUNG FORCEP

SPOUNG FORCEP 10"
FOERSTER SPONGE FORCEPS ARE ESSENTIAL FOR GYNCOLOGIC SURGERY. FEATURING OVAL TIPS WITH SERRATIONS, THE FOERSTER IS USED TO HOLD SPONGES DURING THE PROCEDURE. THE RATCHET HANDLE HELPS TO HOLD THE SPONGES FIRMLY.

KELLEY FORCEP



KELLEY FORCEP

5.5

KELLY FORCEPS THIS INSTRUMENT CAN BE USED AS A CLAMP, HEAT SINK, OR THIRD HAND. HANDLES HAVE TWO-POSITION LOCKING DEVICE. MADE FROM SURGICAL STAINLESS STEEL. SERRATIONS ARE NOT SHARP AND WILL NOT CUT SOFT MATERIALS. SNAP LOCK FOR AN ERGONOMIC HOLD ON SMALL PARTS. CURVED OR STRAIGHT BLADES. 5-1/2" O.L.

ALIS TISSUE FORCEP



ALIS TISSUE FORCEP

REMOVE NECROTIC TISSUE 4.5@


AVAILABLE IN STRAIGHT OR CURVED STYLES
CURVED SCISSORS FOLLOW TISSUE CONTOURS IN SMALL AREAS; STRAIGHT SCISSORS CUT ALONG FLAT SURFACES
MORE ABOUT THIS PRODUCT



KOCHER FORCEP

KOCHER FORCEP

KOCHER HEMOSTATIC FORCEPS, 5.5" (14cm), DELICATE STR TIPS



KERRISON FORCEP UP

**KERRISON FORCEP UP/DOWN
200MM SHAFT LENGTH WITH 40 DEGREE UP BITE
WELL BALANCED STAINLESS STEEL CONSTRUCTION
SPECIFICALLY DESIGNED FOR CUTTING AND REMOVING
SMALL PIECES OF BONE DURING SURGERY
ELIMINATES THE RISK OF CROSS CONTAMINATION
BETWEEN PROCEDURES**



KERRISON FORCEP DOWN

**KERRISON FORCEP UP/DOWN
200MM SHAFT LENGTH WITH 40 DEGREE UP BITE
WELL BALANCED STAINLESS STEEL CONSTRUCTION
SPECIFICALLY DESIGNED FOR CUTTING AND REMOVING
SMALL PIECES OF BONE DURING SURGERY
ELIMINATES THE RISK OF CROSS CONTAMINATION
BETWEEN PROCEDURES**



MARYLAND FORCE

MARYLAND FORCE 13"

THIS COMPLETE RANGE OF SINGLE-USE MONO-POLAR INSTRUMENTATION COVERS ALL NECESSARY FUNCTIONS PERFORMED DURING LAPAROSCOPIC SURGERY. AVAILABLE WITH BOTH A SWITCHABLE RATCHETING HANDLE AND A NON-LOCKING HANDLE, THESE MONO-POLAR INSTRUMENTS OFFER CONSISTENT, PRECISE PERFORMANCE.



DENTAL FORCEP

DENTAL FORCEP No 16

**1ST AND 2ND MOLARS - LOWER. AUTOCLAVABLE,
STAINLESS STEEL FORCEPS WITH A CORN ROW HANDLE**



ASIF VAINCE TRADERS

Certificates



ASIF VAINCE TRADERS



Certificate

CE

Application of council Directive 93/42/EEC OF 14 JUNE, 1993 as updated directive 2007/47/EEC for class Medical Devices

This is certify that the products submitted are: MEDICAL DEVICES CLASS I

"Manufacturers & Exporters of Inert Surgical, Dental, Beauty, Manicure, Pedicure, Holloware, veterinary Instruments and All Sorts of Scissors."

Certificate No: SDC-PK-151212

Manufactured By:

ASIF VAINCE TRADERS

Address to which certificate refers:

Urah Chowk ,Ranger Road, Sialkot 51310 Pakistan.

Comply with the applicable requirements of the directive 93/42/EEC as updated directive 2007/47/EEC the technical file of the devices have been assessed according to the procedure of conformity Assessment described in the Module A, Annexure VII

Limitations

The Manufacturer must inform SDC of any substantial changes occurred in the product or process in order to examen whether this certificate remains valid.

Certificate Issue Date: 20th June, 2021

Date of Initial Registration: 20th May, 2021

Certificate Expiry Date: 19th June, 2024

Re-assessment Date: 19th May, 2024

This certificate of Registration is granted subject to the Regulations approved by the Board.



Scheme Manager

Daniel Thomas

Daniel Thomas



info@stanneg.com
www.stanneg.com



ASIF VAINCE TRADERS



Certificate

The Governing Board of

Stanneg (diensten de cert)

hereby grants to

ASIF VAINCE TRADERS

Certificate No: **SDC-PK-151211**

(hereinafter called the Registered Company) the right to be listed in the Directory of Registered Companies in respect of the service listed below. These services shall be offered by the Registered Company at or from only the address given below in accordance with the quality management system in Compliance with requirement of

ISO 13485:2016

Address to which certificate refers:

Urah Chowk ,Ranger Road, Sialkot 51310 Pakistan.

Approved Scope to which certificate refers:

Manufacturers & Exporters of Inert Surgical, Dental, Beauty, Manicure, Pedicure, Holloware, veterinary Instruments and All Sorts of Scissors.

(Further clarification regarding the Scope of this Certificate and the applicability of ISO 13485:2016 required may be obtained by consulting the organization)

Signed for and behalf of the Board

Certificate Issue Date: 20th June, 2021

Date of Initial Registration: 20th May, 2021

Certificate Expiry Date: 19th June, 2024

Re-assessment Date: 19th May, 2024

This certificate of Registration is granted subject to the Regulations approved by the Board.



Scheme Manager

Daniel Thomas

Daniel Thomas



info@stanneg.com
www.stanneg.com



ASIF VAINCE TRADERS



Certificate

The Governing Board of
Stänneg (diensten de cert)
hereby grants to
ASIF VAINCE TRADERS
Certificate No: SDC-PK-151213

(hereinafter called the Registered Company) the right to be listed in the Directory of Registered Companies in respect of the service listed below. These services shall be offered by the Registered Company at or from only the address given below in accordance with the quality management system in Compliance with requirement of
ISO 9001:2015

Address to which certificate refers:
Urah Chowk ,Ranger Road, Sialkot 51310 Pakistan.

Approved Scope to which certificate refers:
"Manufacturers & Exporters of Inert Surgical, Dental, Beauty, Manicure, Pedicure, Holloware, veterinary Instruments and All Sorts of Scissors.
(Further clarification regarding the Scope of this Certificate and the applicability of ISO 13485:2016 required may be obtained by consulting the organization)
Signed for and behalf of the Board



Certificate Issue Date: 20th June, 2021
Date of Initial Registration: 20th May, 2021

Certificate Expiry Date: 19th June, 2024
Re-assessment Date: 19th May, 2024

This certificate of Registration is granted subject to the Regulations approved by the Board.

Scheme Manager



Daniel Thomas



info@stanneg.com
www.stanneg.com



ASIF VAINCE TRADERS

OFFICE OF DISTRICT HEALTH OFFICER SIALKOT

PHONE: 9250064



D.O.NO. 180

DATE: 26-07-2019

FREE SALE CERTIFICATE

To Whom it May Concern

This is hereby certified, according to the present regulation of Government Of Pakistan the products describe below may be freely distributed in Pakistan and can be exported without any regulatory permission, statutory requirements. These products are lawfully sold in Pakistan and legally distributed around the world .

Manufacturing company is in compliance with ISO-9001(2015), ISO-13485(2016), cGMP (FDA-USA) and Medical devices Class 1

This Certificate being used on the specific request of the member firm.

Name of the Products:

- ✓ Dental Instruments
- ✓ Orthodontic Instruments
- ✓ Surgical Instruments
- ✓ Eye Instruments
- ✓ Electro-Medical Instruments

Name of the Manufacturer,

ASIF VAINCE TRADERS

Head Office: Ranger Road Near Kent Housing Sialkot Pakistan

Sub Office : Commissioner Road Opposite Trauma Center Sialkot Pakistan

On and Behalf of





ASIF VAINCE TRADERS



**U.S. FOOD & DRUG
ADMINISTRATION**

June 23, 2020

Mr. Muhammad Asif
CEO
Asif Vaince Traders
Urah Chowk, Ranger Road
Sialkot – 51310 - Pakistan

FEI: 3016894898

&

Mr. Muhammad Ejaz
Auditor
QA International
1st Floor Shahab Center,
Opp S.I.E
Sialkot – 51310 - Pakistan

Dear Mr. Asif and Mr. Ejaz,

This is to acknowledge receipt of the January 22, 2020 letter from Mr. Muhammad Ejaz (Auditor) certifying the compliance of Asif Vaince Traders (Firm) with the United States Food and Drug Administration (FDA) Quality System Regulation of 1997, which includes the current good manufacturing practice (cGMP) requirements. The Quality System Regulation is set forth in Title 21, Code of Federal Regulations (CFR), Part 820.

The quality system audit report provided by the Auditor states that Asif Vaince Traders manufactures surgical instruments and that a quality system audit was performed on January 13, 2020. The inspection found deficiencies and the Auditor states that a corrective action plan was implemented on January 21, 2020. The Auditor recommends that the Firm be added to the Green List of Import Alert 76-01, "Exempt from Detention Without Physical Examination of Medical Instruments from Pakistan."

FDA has reviewed the audit report of Asif Vaince Traders, including the Quality System Manual, Test Data and the Corrective Action Plan submitted.

Based on our review, Asif Vaince Traders will be placed on the Green List of Import Alert 76-01. The Firm may begin exporting devices to the United States that were manufactured after the consultant certified the Firm's compliance with the cGMP's; however, the Firm's shipments are subject to the guidance outlined in the revised Import Alert 76-01.



Mr. Asif and Mr. Ejaz
Page 2

The FDA may periodically detain and sample devices from the Firm for verification of conformance to the Quality System Regulation. Failure of the sample will result in the Firm being removed from the Green List of Import Alert 76-01 until the Firm is re-inspected and documentation is submitted to the FDA to show compliance with the Quality System Regulation.

The Firm's placement on the Green List of Import Alert 76-01 is limited to medical instrument devices manufactured under the name of Asif Vaince Traders, Urah Chowk, Ranger Road—51310 - Pakistan. In the event the Firm's name and/or manufacturing facility address change, FDA requests that notification be immediately forwarded to this office. A change in the name and/or address of the Firm's manufacturing facility without notifying FDA will result in a re-evaluation of the compliance status of the Firm.

The decision based on the Auditor's certification will remain in effect until such time that FDA is able to visit the Firm for an inspection of the manufacturing facility. During this inspection all corrections and procedures will be evaluated and confirmed. Any new cGMP deviations, or any uncorrected deviations that were previously certified to, may result in a re-evaluation of the compliance status of the Firm, including the possibility of removal from the Green List of Import Alert 76-01. You will be advised of the timing of FDA's inspection schedule.

If the Firm has not conducted a Quality System audit in the past two years, we request that a Quality System audit be conducted within 6 months of receiving this letter. A copy of the Firm's most recent audit should be submitted to FDA for review. The Firm has a responsibility to conduct periodic Quality System audits to ensure conformance with the Quality System regulation.

The audit report should address, at a minimum, the applicable elements of the Quality System Regulation including the following information, as appropriate. This should not be considered an all-inclusive list and additional information may be included.

- | | |
|-----------------------------------|--------------------------------|
| o Current Audit Summary | o Facilities |
| o Audit Follow-up Recommendations | o Supplier Control |
| o Corrective Action Plan | o Specifications |
| o Device Master Record | o Production Equipment |
| o Device History Record | o Cleaning and Sanitation |
| o Material Analysis | o Personal Hygiene |
| o Material Hardness Tests | o Training |
| o Material Corrosion Tests | o Hazardous Materials Handling |
| o Quality Manual | o Receiving, Storage, Shipping |
| o Calibration | o Traceability and Recall |
| o Internal Audits | o Consumer Complaints/MDRs |
| o External Audits | o Pest Control |

All manufacturers exporting surgical instruments to the United States should use stainless steel meeting the latest version of the Standard Specification for Wrought Stainless Steels for Surgical Instruments, ASTM standard F-899-12B. Please assure that the Firm's documents and



Mr. Asif and Mr. Ejaz
Page 3

requirements conform to ASTM standard F-899-12B.

Establishments that are involved in the production and distribution of medical devices intended for use in the United States are required to register and list the devices annually with the FDA. This registration and listing process may be completed electronically. For more information and to complete the process please go to:

<http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/RegistrationandListing/default.htm>.

Electronic submission documents should be emailed to **cdrhocpakistanaudit@fda.hhs.gov**.

Cover letters contained in the electronic submission should be addressed to:

Ms. Deniz B. Mackey, Assistant Director
Imports and Registration & Listing Team
U.S. Food and Drug Administration - CDRH
Office of Regulatory Programs
Division of Regulatory Programs 2
White Oak Building 66
10903 New Hampshire Avenue
Silver Spring, Maryland 20993 USA

Please reference your Facility Establishment Number (FEI), 3016894898, in future correspondence and in the registration process.

If you have any questions regarding this correspondence, or need further assistance, please contact Raffi Papazian at raffi.papazian@fda.hhs.gov or (240) 402-1670.

Sincerely yours,

 Digitally signed by
Tanisha L. Hithé-S
Date: 2020.06.23
20:04:55 -04'00'

Deniz B. Mackey
Assistant Director
Imports and Registration & Listing Team
Division of Regulatory Programs 2
Office of Regulatory Programs
Center for Devices and Radiological Health